IVD In vitro Diagnostic Use Only

PR-CVDCAgS1

REF PR-CVDCAgS5

PR-CVDCAgS20

Σ 1, 5 or 20

Please Read the Instruction Sheet Carefully Before Performing the Test.

INTENDED USE

ProDetect[®] COVID-19 Antigen Rapid Self-Test (Saliva) is a rapid lateral flow chromatographic immunoassay for the qualitative detection of SARS-CoV-2 nucleocapsid protein antigens in human saliva specimens. The test serves as an aid in the diagnosis of coronavirus infection disease (COVID-19) for individuals suspected of SARS-CoV-2 infection. ProDetect[®] COVID-19 Antigen Rapid Self-Test (Saliva) is used in conjunction with clinical presentation and the results of other laboratory tests. It is intended for home use only.

SUMMARY AND EXPLANATION OF THE TEST

Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) was identified as the causative agent for the outbreak of pneumonia cases in early January 2020. It is associated with common symptoms including fever, cough, fatigue, shortness of breath or breathing difficulties, and loss of smell and taste. While most people have mild symptoms, some people develop acute respiratory distress syndrome (ARDS). It starts with mild symptoms for about a week and then progresses to rapid deterioration and ARDS requiring advanced life support.

ProDetect[®] COVID-19 Antigen Rapid Self-Test (Saliva) detects SARS-CoV-2 nucleocapsid protein antigens. During the acute phase of infection, an antigen is commonly detectable in the upper respiratory specimens.

PRINCIPLE OF THE TEST

ProDetect[®] COVID-19 Antigen Rapid Self-Test (Saliva) is a qualitative membrane-based immunoassay for the detection of SARS-CoV-2 nucleocapsid protein antigens in human saliva specimen. During the testing, the specimen reacts with SARS-CoV-2 nucleocapsid protein antibody-coated particles in the test line region of the test device. The mixture then migrates upward on the membrane chromatographically by capillary action and reacts with SARS-CoV-2 nucleocapsid protein antibody in test line region. A colored line will appear in test line region if specimen contains SARS-CoV-2 antigens. If the specimen does not contain antigens to SARS-CoV-2, no colored line will appear in the test line region, indicating a negative result. To serve as a procedure control, a colored line will appear at the Control Region (C), indicating that the test has been performed properly.

REAGENTS AND MATERIALS SUPPLIED

Catalog Number	Item
PR-CVDCAgS1/ PR-CVDCAgS5/ PR-CVDCAgS20	One box consists of the following components: 1. 1/5/20 Individually sealed pouches, each containing: <ul style="list-style-type: none"> One test device One desiccant pouch 2. 1/5/20 collection device(s), each consisting of: <ul style="list-style-type: none"> One funnel One collection tube 3. 1/5/20 biosafety bag(s) 4. 1/5/20 single use buffer(s) 5. 1 leaflet with instruction for use

WARNINGS AND PRECAUTIONS

- This test device is for *in vitro* diagnostic use only.
- This instruction for use must be read completely before performing the test. Failure to follow directions in the instruction for use may yield inaccurate test results.
- Do not eat, drink, or smoke in the area where the specimens or test devices are handled.
- Clean up spills thoroughly with appropriate disinfectants.
- Avoid splashing or aerosol formation while collecting the sample.
- Discard after first use. The test device cannot be used more than once.
- Do not use the test device beyond the expiry date.
- Do not use the test device if the pouch is punctured or not well sealed.
- Do not mix and interchange different specimens.
- Keep out of reach of children.
- Do not read after 20 minutes.
- Do not swallow.
- Discard of all specimens and used test device into biosafety bag.

STORAGE

- The test device is stable until the indicated expiry date when stored in the sealed pouch at 2-30°C and protected from direct sunlight, moisture, and heat.
- DO NOT FREEZE.
- The test device must remain in the sealed pouch until use.

SAMPLE COLLECTION AND PREPARATION

Specimen Collection and Extraction

Saliva specimen should be collected using the collection device provided with the test device. Do not use other collection devices. Saliva collected at any time of the day may be used. Specimens should be tested immediately after collection if possible.

Prior to collection of saliva, you must not place anything into your mouth including food, drink, gum, or tobacco products for at least 10 minutes before collection.

Step 1:

- Attach the funnel to the collection tube prior to sample collection.
- Deeply cough 3-5 times and spit saliva into the collection device.
Note: It is recommended to collect the first saliva after deep coughing in the morning.
- Repeat the steps above for sample collection until the saliva reaches the scale line indicated on the collection tube.

Step 2:

- Remove the funnel from collection tube.
- Empty the entire content of the disposable buffer into the collection tube.
- Place the cap securely on the collection tube. Gently squeeze the collection tube for 10-15 times to mix the sample with the buffer thoroughly.




TEST PROCEDURE


For video demonstration, please scan the QR code below:

Bring the test device, buffer, and sample to room temperature (15-30°C) prior to testing. Best results will be obtained if the test is performed immediately after opening the foil pouch.

- Remove device from the foil pouch (use within 1 hour) and place it on a flat and clean surface.
- Transfer extracted sample to the test device.
- Wait for the colored line(s) to appear. Read the results at 15 minutes.
- Do not read results after 20 minutes.
- Discard the used test device into the biosafety bag and dispose the biosafety bag.

INTERPRETATION AND REPORTING OF RESULTS

Positive (+)

Two colored lines appear at test line (T) and control line (C). It indicates a positive result for the SARS-CoV-2 antigens in the specimen.
Note: Color intensity of the line appearing in the test (T) region may vary depending on the SARS-CoV-2 antigen level in the specimen. Therefore, any shade of color in the test (T) region is to be considered as a positive result.

Procedure in reporting positive result via MySejahtera:

- If positive result is obtained, notify Ministry of Health by either using self-notification at *MySejahtera* mobile application or notify the District Health Office.
- Covid-19 positive individual should contact the nearest COVID-19 Assessment Center (CAC) (<http://covid-19.moh.gov.my/hotline>) to inform and get advice for further action to be taken.

1. Scan the QR code above to link you to MySejahtera Helpdesk page.
2. Select whether you are reporting for yourself or for your dependent.
3. Once select, you are required to update individual particulars (i.e name, mobile, identification number, and email address) as required.
4. Stay on the Helpdesk page and continue to update with your home address.
5. Stay on the Helpdesk page and select one of the options ("hospital/clinic"; "online"; "pharmacy"; "corporate", or "others") as to where you've obtained the self-test kits from.
6. Stay on the Helpdesk page and select one of the options ("saliva"; "nasal", or "others") as to what was tested with the self-test kits.
7. To submit the self-test results, on the following page, select one of the options ("positive"; "negative", or "invalid").

Negative (-)

The colored line appears at the control line (C) only. It indicates that the concentration of the SARS-CoV-2 antigen is zero or below the detection limit of the test.

Invalid

Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test device immediately and contact your local distributor.

LIMITATIONS OF THE TEST

1. The test procedure must be carefully followed when testing the presence of SARS-CoV-2 nucleocapsid protein antigens in human saliva specimen from suspected individuals. Failure to follow the instructions may cause incorrect results.
2. The performance of ProDetect® COVID-19 Antigen Rapid Self-Test (Saliva) was validated using the procedures provided in this instruction for use. Modification to these procedures may cause inaccurate performance.
3. ProDetect® COVID-19 Antigen Rapid Self-Test (Saliva) is for *in vitro* diagnostic use only. It is used for the detection of SARS-CoV-2 nucleocapsid protein antigens in human saliva specimens as an aid in the diagnosis of individuals suspected of SARS-CoV-2 infection in conjunction with clinical symptoms and the results of other laboratory tests. This qualitative test does not provide quantitative value and cannot determine the rate of increase in the concentration of SARS-Cov-2 nucleocapsid protein antigens.
4. This test will only indicate the presence of SARS-CoV-2 antigens in the specimen and should not be used as the sole criteria for the diagnosis of SARS-CoV-2 infection. The results obtained with the test should be considered with other clinical findings from other laboratory tests and evaluations.
5. If the test result is negative or non-reactive but clinical symptoms persist, it is recommended that you test again or test with a molecular diagnostic device to rule out infection.
6. Negative results may be obtained if the titre of the novel coronavirus antigens in the sample is lower than the minimum detection limit of the test. However, they do not rule out SARS-CoV-2 infection, particularly in

individuals who have been in contact with the virus. Follow-up testing with a molecular diagnostic should be considered to rule out infection.
 7. Infection with non-SARS-CoV-2 coronavirus strains or other interference factors may show in positive results with ProDetect® COVID-19 Antigen Rapid Self-Test (Saliva).

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity:

ProDetect® COVID-19 Antigen Rapid Self-Test (Saliva) was compared with a leading commercial PCR test device. The results showed that the ProDetect® COVID-19 has high sensitivity and specificity.

ProDetect® COVID-19 Antigen Rapid Self-Test (Saliva)	Method	PCR		Total Results
	Results	Positive	Negative	
Antigen Rapid Self-Test (Saliva)	Positive	91	2	93
	Negative	10	303	313
Total Result		101	305	406
Relative Sensitivity		90.1% (95%CI: 82.5% to 95.1%)		
Relative Specificity		99.3% (95%CI: 97.7% to 99.9%)		
Accuracy		97.0% (95%CI: 94.9% to 98.5%)		

CI: Confidence Intervals

Specificity with other viral strains

ProDetect® COVID-19 Antigen Rapid Self-Test (Saliva) has been tested with other viral strains. The results showed no discernible line at the test-line regions at the following concentrations:

Virus Type	Concentration (TCID ₅₀ /ml)
Adenovirus type 3	3.16 x 10 ⁴
Adenovirus type 7	1.58 x 10 ⁵
Human coronavirus OC43	1.00 x 10 ⁶
Human coronavirus 229E	5.00 x 10 ⁵
Human coronavirus HKU1	1.00 x 10 ⁶
Human coronavirus NL63	1.00 x 10 ⁶
Influenza A H1N1	3.16 x 10 ⁵
Influenza A H3N2	1.00 x 10 ⁵
Influenza B	3.16 x 10 ⁶
Measles	1.58 x 10 ⁴
Mumps	1.58 x 10 ⁴
Parainfluenza virus 2	1.58 x 10 ⁷
Parainfluenza virus 3	1.58 x 10 ⁸
Respiratory syncytial virus	8.89 x 10 ⁴

Precision (Intra-assay & Inter-assay)

Precision was determined by using three specimens of COVID-19 standard control. Three different lots of ProDetect® COVID-19 Antigen Rapid Self-Test (Saliva) have been tested using negative, weak SARS-CoV-2 antigen and strong SARS-CoV-2 antigen. Ten replicates of each level were tested each day for 3 consecutive days. The specimens were correctly identified >99% of the time.

CROSS REACTIVITY

The following organisms were tested at the concentration of 1.0 x 10⁸ org/ml using the test device. All were found to be negative.

<i>Arcanobacterium</i>	<i>Pseudomonas aeruginosa</i>
<i>Candida albicans</i>	<i>Staphylococcus aureus sub aureus</i>
<i>Corynebacterium</i>	<i>Staphylococcus epidermidis</i>
<i>Escherichia coli</i>	<i>Streptococcus pneumoniae</i>
<i>Moraxella catarrhalis</i>	<i>Streptococcus pyogenes</i>
<i>Neisseria lactamica</i>	<i>Streptococcus salivarius</i>
<i>Neisseria subflava</i>	<i>Streptococcus sp. Group F</i>

Interfering Substances


The following substances were tested with ProDetect® COVID-19 Antigen Rapid Self-Test (Saliva) and no interference was observed:

Caffeine	1mg/ml
Dexamethasone	0.8mg/ml
Flunisolide	6.8ng/ml
Milk	11.2%
Mouthwash	2.0%
Mucin	50µg/ml
Mupirocin	12mg/ml
Orange Juice	100.0%
Oxymetazoline	0.6mg/ml
Phenylephrine	12.0 g/ml
Rebetol	4.5µg/ml
Relenza	282ng/ml
Tamiflu	1.1µg/ml
Tea	33.3mg/ml
Tobramycin	2.43mg/ml

REFERENCES

1. Huang, C. et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. *Lancet* 395, 497–506 (2020).
2. Steph H Tan, Orchid Allcock, Mari Armstrong-Hough, Anne L Wyllie. Saliva as a gold-standard sample for SARS-CoV-2 detection. *The Lancet Respiratory Medicine*, 562-564 (2021).
3. *Emerging respiratory viruses, including COVID-19: methods for detection, prevention, response and control.* (2020). World Health Organization. <https://openwho.org/courses/introduction-to-ncov>

INDEX OF SYMBOLS USED

 Consult instructions for use	 Do not re-use
 Catalog number	 No. of tests
 Batch code	 Temperature limitation
 Date of manufacture	 Use by
 Manufacturer	

 Medical Innovation Ventures Sdn. Bhd. (Co. No. 988633-U),
1st Floor, Plot 88f, Lintang Bayan Lepas 10,
Bayan Lepas Industrial Park, Phase 4,
11900 Bayan Lepas, Pulau Pinang, Malaysia.
Phone: +604-3052730 Fax: +604-3052730
Email: sales@mediven.com.co

 and  logos are trademarks of

Medical Innovation Ventures Sdn. Bhd.

PRCVDCAgS- Rev. 0 English Version
Effective: 2021- Nov Issued: 2021 Nov

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RUJUKAN

- Huang, C. et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. *Lancet* 395, 497–506 (2020).
- Steph H Tan, Orchid Allicock, Mari Armstrong-Hough, Anne L Wyllie. Saliva as a gold-standard sample for SARS-CoV-2 detection. *The Lancet Respiratory Medicine*, 562-564 (2021).
- Emerging respiratory viruses, including COVID-19: methods for detection, prevention, response and control. (2020). World Health Organization. <https://openwho.org/courses/introduction-to-ncov>

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