TEPAT COVID-19 Ag Rapid Self-Test (Saliva)

IVD

In vitro Diagnostic Use Only

REF

TP-CVDAgSL

Σ 1,5 or 25

Please Read the Instruction Sheet Carefully Before Performing the Test.

INTENDED USE

TEPAT COVID-19 Ag Rapid Self-Test (Saliva) is a rapid lateral flow chromatographic immunoassay for the qualitative detection of SARS-CoV-2 nucleocapsid protein antigens in human saliva specimens. The test serves as an aid in the diagnosis of COVID-19 infection for suspected individuals. TEPAT COVID-19 Ag Rapid Self-Test (Saliva) is used in conjunction with clinical presentation and the results of other laboratory tests. It is intended for home use only.

SUMMARY AND EXPLANATION OF THE TEST

Previously known as the 2019-novel coronavirus (2019-nCoV), the novel virus is now known as the severe acute respiratory syndrome coronavirus-2 (SARS-CoV-2) and has been identified as the primary agent of an outbreak of pneumonia cases. In early 2020, people all over the world are alarmed as the emergence of this novel coronavirus has posed significant public health and governance challenges. SARS-CoV-2 is an enveloped, positive-sense, single-stranded RNA *Betacoronavirus* closely related to the viruses that caused the 2002–2004 SARS epidemics and the 2012 Middle East respiratory disease (MERS) outbreaks.

The SARS-CoV-2 infection is associated with common symptoms including fever, cough, fatigue, shortness of breath or breathing difficulties, and loss of smell and taste. While most people have mild symptoms, some people develop acute respiratory distress syndrome (ARDS). It begins with mild symptoms for about a week and then progresses to rapid deterioration and ARDS, necessitating extensive life support.

TEPAT COVID-19 Ag Rapid Self-Test (Saliva) detects SARS-CoV-2 nucleocapsid protein antigens. During the acute phase of infection, antigen testing is commonly used on specimens collected within 7 days of the onset of symptoms in a suspected population.

PRINCIPLE OF THE TEST

TEPAT COVID-19 Ag Rapid Self-Test (Saliva) is a qualitative membrane-based immunoassay designed to detect SARS-CoV-2 nucleocapsid protein antigens in human saliva specimens. The specimen interacts with antibody-coated SARS-CoV-2 nucleocapsid protein particles in the test (T) line region of the test device during the testing. The mixture then migrates upward on the membrane chromatographically by capillary action and interacts with SARS-CoV-2 nucleocapsid protein antibody in T line region. If the specimen contains SARS-CoV-2 antigens, a colored line will show in T line region. If no antigens against SARS-CoV-2 in the specimen, no colored line will show in the T line region, indicating a negative result. A colored line will appear at the control (C) line region to serve as a procedure control, indicating that the test has been performed properly.

REAGENTS AND MATERIALS SUPPLIED

- 1/5/25 individually sealed pouches, each containing:
 - Test device
 - Desiccant pouch
- 1/5/25 collection device(s), each consisting of
 - Collection funnel
 - Extraction tube
- 1/5/25 Extraction buffer(s)
- 1/5/25 Biosafety bag(s)
- Work station (only available for 5/25 tests/box)
- 1 leaflet with instruction for use

MATERIALS NEEDED BUT NOT SUPPLIED

1. Timer
2. Disinfectants

WARNINGS AND PRECAUTIONS

1. This kit is for *in vitro* diagnostic use only.
2. This instruction for use must be read completely before performing the test. Failure to follow directions in the instruction for use may yield inaccurate test results.
3. Bloody specimens must not be used for testing.
4. Do not touch the membrane and specimen well.
5. Clean up spills thoroughly with appropriate disinfectants.
6. Discard after first use. The test device cannot be used more than once.
7. Do not use the test device beyond the expiry date.
8. Do not use the test device if the pouch is punctured or not well sealed.
9. Do not mix and interchange different specimens.
10. Discard of all specimens and used test device into biosafety bag.
11. Do not eat, drink or smoke in the area where the specimens or kits are handled.
12. All specimens should be considered potentially hazardous and handled in the same manner as an infection agent.

STORAGE

1. The test device is stable until the indicated expiry date when stored in the sealed pouch at 2-30°C and protected from direct sunlight, moisture, and heat.
2. DO NOT FREEZE.
3. The test device must remain in the sealed pouch until use.

PURCHASE RECORD AND RESULTS REPORTING

1. Scan the QR code above to download MyDocLab app.
2. Create an account and edit your profile.
3. Register the purchased self-test kit device by selecting "Purchase Record" and fill in the required information for submission.
4. For self-test result submission, select "Self-Test".

- Step 1: Validate your purchased self-test kit by choosing one of the available options.
- Step 2: Once your self-test kit is verified as genuine, click "Next" to begin testing.
- Step 3: Select either "Individual" or "Corporate" as your submission type. Then, click "Next" to proceed. Concerning "Corporate," you may want to scan Corporate QR code/input the Corporate code provided by your corporate management upon clicking "Next".
- Step 4: You may begin the countdown timer by clicking "Start" once the specimen has been added to the self-test kit.
- Step 5: To submit your self-test results, select one of the options ("positive"; "negative", or "invalid") and upload a snapshot of your result.

5. To report your self-test result on MySejahtera app, go to the "Menu" and click the previously reported result. Following that, select "Submission to MySejahtera" to be directed to MySejahtera Helpdesk page.

- Step 1: Select if you are reporting for yourself or a dependent.
- Step 2: Once selected, you are required to update your individual particulars (i.e name, mobile, identification number, and email address) as required.
- Step 3: Continue to update your home address on the Helpdesk page.
- Step 4: Remain on the Helpdesk page and select one of the options ("hospital/clinic"; "online"; "pharmacy"; "corporate", or "others") as to where you've obtained the self-test kits from.
- Step 5: Remain on the Helpdesk page and select one of the options ("saliva"; "nasal", or "others") for the type of test performed using the self-test kits.
- Step 6: On the next page, pick one of the choices for submitting the self-test results ("positive"; "negative", or "invalid").

SPECIMEN COLLECTION AND PREPARATION

Bring the test device and components of the kit to room temperature (15-30°C) prior to testing. The best results will be obtained if the test is performed immediately after opening the foil pouch.

1. Prior to testing, wash hands with soap and water thoroughly.
2. Do not place anything into mouth including food, drink, gum, or tobacco products for at least 30 minutes before saliva collection.
3. Avoid contact to infections like severe mouth ulcers and bronchitis as this may affect saliva collection.
4. Rinse mouth with running water 30 minutes before saliva collection.

Specimen Collection

1. Attach the collection funnel to the extraction tube prior to specimen collection.
2. Deeply cough 3-5 times and spit saliva into the extraction tube until the specimen could reach the scale 300µL on the wall of extraction tube excluding bubbles is sufficient.
3. Repeat the steps (1-2) for specimen collection until the saliva reaches the scale line indicated on the extraction tube.

Note: For best results, saliva specimen collected in the morning before mouth rising, eating or drinking is recommended.

4. When the specimen is seen reaching the scale line indicated on the extraction tube, remove the collection funnel and place the extraction tube in the work station.

Specimen Processing

1. Empty the entire content of the extraction buffer into the extracted tube containing saliva specimen.
2. Place the cap securely on the extraction tube and gently squeeze the extraction tube for 10-15 times to mix the specimen with the extraction buffer thoroughly.

Note: The processed specimen should be tested within 1 hour.


TEST PROCEDURE

For video demonstration, please scan the QR code below:



1. Remove test device from the foil pouch and place it on a flat and clean surface.
2. Add 2 drops of processed specimen to the specimen well of the test device.
3. Wait for the colored line(s) to appear. Read the results at 20 minutes.
4. Do not read results after 30 minutes.
5. Discard the used test device into the biosafety bag and dispose the biosafety bag.

Specimen Collection

Specimen Processing

INTERPRETATION AND REPORTING OF RESULTS

Positive (+)

Two colored lines appear at test (T) line and control (C) line. It indicates a positive result for the SARS-CoV-2 antigens in the specimen.

Note: Color intensity of the line appearing in the test (T) region may vary depending on the SARS-CoV-2 antigen level in the specimen. Therefore, any shade of color in the test (T) region is to be considered as a positive result.


- If positive result is obtained, notify Ministry of Health by either using self-notification at MySejahtera mobile application or notify the District Health Office.
- Covid-19 positive individual should contact the nearest COVID-19 Assessment Center (CAC) (<http://covid-19.moh.gov.my/hotline>) to inform and get advice for further action to be taken.

Negative (-)

The colored line appears at the control (C) line only. It indicates that the concentration of the SARS-CoV-2 antigen is zero or below the detection limit of the test.

Invalid

Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test device. If the problem persists, discontinue using the test device immediately and contact your local distributor.

LIMITATIONS OF THE TEST

- The test procedure must be carefully followed when testing the presence of SARS-CoV-2 nucleocapsid protein antigens in human saliva specimen from suspected individuals. Failure to follow the instructions may cause incorrect results.
- The performance of TEPAT COVID-19 Ag Rapid Self-Test (Saliva) was validated using the procedures provided in this instruction for use. Modification to these procedures may cause inaccurate performance.
- TEPAT COVID-19 Ag Rapid Self-Test (Saliva) is for *in vitro* diagnostic use only. It is used for the detection of SARS-CoV-2 nucleocapsid protein antigens in human saliva specimens as an aid in the diagnosis of individuals suspected of SARS-CoV-2 infection in conjunction with clinical symptoms and the results of other laboratory tests. This qualitative test does not provide quantitative value and cannot determine the rate of increase in the concentration of SARS-CoV-2 nucleocapsid protein antigens.
- This test will only indicate the presence of SARS-CoV-2 antigens in the specimen and should not be used as the sole criteria for the diagnosis of SARS-CoV-2 infection. The results obtained with the test should be considered with other clinical findings from other laboratory tests and evaluations.
- If the test result is negative or non-reactive but clinical symptoms persist, it is recommended that you test again or test with a molecular diagnostic device to rule out infection.
- Negative results may be obtained if the titre of the novel coronavirus antigens in the specimen is lower than the minimum detection limit of the test. However, they do not rule out SARS-CoV-2 infection, particularly in individuals who have been in contact with the virus as the expression of the antigen might be low during the early stages of infection. To rule out infection, follow-up testing using a molecular diagnosis should be considered.
- Infection with non-SARS-CoV-2 coronavirus strains or other interference factors may show in positive results with TEPAT COVID-19 Ag Rapid Self-Test (Saliva).
- Blood or excess mucus inside the saliva specimen might cause interference and false positive results.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity

TEPAT COVID-19 Ag Rapid Self-Test (Saliva) was compared with a leading commercial PCR test device. The results showed that the TEPAT COVID-19 Ag Rapid Self-Test (Saliva) has high sensitivity and specificity.

TEPAT COVID-19 Ag Rapid Self-Test (Saliva)	Method	PCR		Total Results
	Result	Positive	Negative	
	Positive	96	3	
	Negative	4	297	301
Total Result		100	300	400
Relative Sensitivity		96.00% (95%CI: 90.07% to 98.90%)		
Relative Specificity		99.00% (95%CI: 97.11% to 99.79%)		
Accuracy		98.25% (95%CI: 96.43% to 99.29%)		

CI: Confidence Intervals

Precision (Intra-assay & Inter-assay)

Precision was determined by using three specimens of COVID-19 standard control. Three different lots of TEPAT COVID-19 Ag Rapid Self-Test (Saliva) have been tested using negative, weak SARS-CoV-2 antigen and strong SARS-CoV-2 antigen. Five replicates of each level were tested each day for 3 consecutive days. The specimens were correctly identified >99% of the time.

CROSS REACTIVITY

The following pathogens were tested at concentration 10^5 TCID₅₀/ml and 10^6 TCID₅₀/ml, respectively using the test device. All pathogens showed no cross reactivity with the test device.

Conc.	List of Pathogens Tested	
10^5 TCID ₅₀ /ml	Adenovirus type 3 Epstein-Barr virus Enterovirus CA16 HCOV-HKU1 Human metapneumovirus Human coronavirus 229E Human coronavirus NL63 Human coronavirus OC43 H7N9 H5N1	Influenza A H1N1 2009 Influenza A H3N2 Influenza B victoria strain Influenza B Y strain MERS coronavirus Measles virus Mumps virus Parainfluenza virus, type2 Rhinovirus Respiratory syncytial virus
10^6 TCID ₅₀ /ml	<i>Bordetella parakeratosis</i> <i>Bordetella pertussis</i> <i>Candida albicans</i> <i>Chlamydia pneumoniae</i> <i>Group A streptococci</i> <i>Legionella pneumophila</i> <i>Mycobacterium tuberculosis</i> <i>Mycoplasma pneumoniae</i> <i>Pneumocystis jirovecii</i> <i>Staphylococcus aureus</i> <i>Streptococcus pneumoniae</i>	

Interfering Substances

The following substances were tested with TEPAT COVID-19 Ag Rapid Self-Test (Saliva) and no interference was observed:

Substance	Conc.	Substance	Conc.
Afrin	15%	Menthol	15%
Compound Benzoin Gel	1.5mg/ml	Mucin	0.5%
Cromolyn Glycate	15%	Mupirocin	10mg/ml
Chloramphenicol	3ug/ml	Naphazoline Hydrochloride Nasal Drops	15%
Deoxyepinephrine Hydrochloride	15%	Oseltamivir	5mg/ml
Erythromycin	3ug/ml	Tetracycline	3ug/ml
Fluticasone Propionate Spray	15%	Tobramycin	5%
Ibuprofen	1mg/ml	Whole Blood	4%

REFERENCES

- Zheng J. (2020). SARS-CoV-2: an Emerging Coronavirus that Causes a Global Threat. *International journal of biological sciences*, 16(10), 1678–1685. <https://doi.org/10.7150/ijbs.45053>.
- Riva, L., Yuan, S., Yin, X. *et al.* Discovery of SARS-CoV-2 antiviral drugs through large-scale compound repurposing. *Nature* 586, 113–119 (2020). <https://doi.org/10.1038/s41586-020-2577-1>
- Steph H Tan, Orchid Allicock, Mari Armstrong-Hough, Anne L Wyllie. Saliva as a gold-standard sample for SARS-CoV-2 detection. *The Lancet Respiratory Medicine*, 562-564 (2021).
- Emerging respiratory viruses, including COVID-19: methods for detection, prevention, response and control.* (2020). World Health Organization. <https://openwho.org/courses/introduction-to-ncov>

- Consult instructions for use
- Manufacturer
- For *in vitro* diagnostic use only
- No. of tests
- Catalog number
- Use by
- Batch code
- Temperature limitation
- Date of manufacture
- Do not re-use

Medical Innovation Ventures Sdn. Bhd. (Co. No. 988633-U),
1st Floor, Plot 88f, Lintang Bayan Lepas 10,
Bayan Lepas Industrial Park, Phase 4,
11900 Bayan Lepas, Penang, Malaysia.
Phone: +604-305 2730 Fax: +604-305 2730
Email: sales@mediven.com.co

MEDIVEN
MEDICAL INNOVATION VENTURES

and

TEPAT

logos are trademarks of

Medical Innovation Ventures Sdn. Bhd.

PRCVDAGSL - Rev. 0 English Version
Effective: 2022-Feb Issued: 2022-Feb

TEPAT COVID-19 Ag Rapid Self-Test (Saliva)

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Kepekatan	96.00% (95%CI: 90.07% to 98.90%)
Kekhususan	99.00% (95%CI: 97.11% to 99.79%)
Ketepatan	98.25% (95%CI: 96.43% to 99.29.%)

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RUJUKAN

- Zheng J. (2020). SARS-CoV-2: an Emerging Coronavirus that Causes a Global Threat. *International journal of biological sciences*, 16(10), 1678–1685. <https://doi.org/10.7150/ijbs.45053>.
- Riva, L., Yuan, S., Yin, X. *et al.* Discovery of SARS-CoV-2 antiviral drugs through large-scale compound repurposing. *Nature* 586, 113–119 (2020). <https://doi.org/10.1038/s41586-020-2577-1>
- Steph H Tan, Orchid Alicock, Mari Armstrong-Hough, Anne L Wyllie. Saliva as a gold-standard sample for SARS-CoV-2 detection. *The Lancet Respiratory Medicine*, 562-564 (2021).
- Emerging respiratory viruses, including COVID-19: methods for detection, prevention, response and control.* (2020). World Health Organization. <https://openwho.org/courses/introduction-to-ncov>

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