SARS-CoV-2 Antigen Saliva Lolly Test Package Insert

| | | | |
|----------------|----------|-----------------------|---------|
| REF | COVG-603 | Sample: | Saliva |
| Version | Z | Effective Date | 2022.03 |

INTENDED USE

The SARS-CoV-2 Antigen Saliva Lolly Test is used for in vitro qualitative determination of SARS-CoV-2 antigens in human saliva samples. It can be used for rapid investigation of suspected COVID-19 cases and can be used as a reconfirmation method for nucleic acid detection in discharged cases. A positive test result indicates that the sample contains SARS-CoV-2 antigen. A negative test result does not rule out the possibility of infection.

The SARS-CoV-2 Antigen Saliva Lolly Test is intended for home use laymen in a non-traditional setting (such as offices, sporting events, airports, schools, etc.). The test result of this kit is for clinical reference only. It is recommended to conduct a comprehensive analysis of the condition based on the patient's clinical manifestations and other laboratory tests. Antigen testing is typically used in the acute phase of infection, when samples are tested within seven days of the onset of symptoms in a suspected population.

SUMMARY

The novel coronaviruses belong to the beta genus. COVID-19 is an acute infectious disease of the respiratory tract. Currently, patients infected with the novel coronavirus are the main source of infection. Infected people without symptoms can also infect others. The incubation period is 1 to 14 days, usually 3 to 7 days. The main symptoms are fever, fatigue and a dry cough. Nasal congestion, runny nose, sore throat, muscle pain, and diarrhea occur in some cases.

PRINCIPLE

The SARS-CoV-2 Antigen Saliva Lolly Test is a qualitative, lateral flow immunoassay for the detection of the N protein of SARS-CoV-2 in human saliva. In this test, antibody specific to the N protein of SARS-CoV-2 is separately coated on the test line region of the test midstream. During testing, the extracted specimen reacts with the antibody to N protein of SARS-CoV-2 that are coated onto particles. The mixture migrates up to the membrane to react with the antibody to N protein of SARS-CoV-2 on the membrane and generate one colored line in the test region. The presence of this color line of test region indicates a positive result. To serve as a procedural control, a colored line will always appear in the control region if the test has performed properly.

PRECAUTIONS

- Please read all the information in this package insert before performing the test.
- For self-testing in vitro diagnostic use only. Do not use after expiration date.
- The test midstream should remain in the sealed pouch until ready to use.
- All samples should be considered potentially hazardous and handled in the same manner as an infection agent.
- Avoid using bloody samples.
- Test for children and young people should be tested with an adult supervision.
- Discard the test midstream into the biohazard waste bag and dispose it according to local regulations.
- Do not reuse.
- Wash hands thoroughly before and after handling.
- Please ensure that an appropriate amount of samples are used for testing. Too much or too little sample size may lead to deviation of results.
- 30 minutes before saliva collection procedure, DO NOT eat, drink, smoke or chew gum.
- 30 minutes before saliva collection procedure, use running water to clean mouth for a few times.

STORAGE AND STABILITY

- Store as packaged in the sealed pouch at room temperature or refrigerated (2-30 °C).
- The test is stable through the expiration date printed on the sealed pouch. Do not use beyond the expiration date.
- The test must remain in the sealed pouch until use.
- **DO NOT FREEZE.**

PACKAGE SPECIFICATIONS

1 test/pack, 5 tests/pack, 25 tests/pack, 50 tests/pack, 100 tests/pack

MATERIALS PROVIDED

- i) Materials required and provided:
 - SARS-CoV-2 Antigen Saliva Lolly Test midstream
 - Package Insert
 - Qualification Certificate

- Biohazard bag

Note: Components of different batches cannot be mixed.

ii) Materials required but not provided:

The timer and disinfection products, such as hand sanitizer, rubbing alcohol, soap, etc.

PURCHASE RECORD AND RESULTS REPORTING

1. Scan the QR code above to download MyDocLab app.
2. Create an account and edit your profile.
3. Scan the QR code printed on the purchased self-test kit device to register by selecting "Purchase Record".
4. For self-test result submission, select "Self-Test".
 - Step 1: Validate your purchased self-test kit by choosing one of the available options.
 - Step 2: Once your self-test kit is verified as genuine, click "Next" to begin testing.
 - Step 3: Select either "Individual" or "Corporate" as your submission type. Then, click "Next" to proceed. Concerning "Corporate", you may want to scan Corporate QR code/input the Corporate code provided by your corporate management upon clicking "Next".
 - Step 4: You may begin the countdown timer by clicking "Start" once the specimen has adequately absorbed by the self-test kit.
 - Step 5: To submit your self-test results, select one of the options ("positive"; "negative", or "invalid") and upload a snapshot of your result.
5. To report your self-test result on MySejahtera app, go to the "Menu" and click previously reported result. Following that, select "Submission to MySejahtera" to be directed to MySejahtera Helpdesk page.
 - Step 1: Select if you are reporting for yourself or a dependent.
 - Step 2: Once selected, you are required to update your individual particulars (i.e name, mobile, identification number, and email address) as required.
 - Step 3: Continue to update your home address on the Helpdesk page.
 - Step 4: Remain on the Helpdesk page and select one of the options ("hospital/ clinic"; "online"; "pharmacy"; "corporate", or "others") as to where you've obtained the self-test kits from.
 - Step 5: Remain on the Helpdesk page and select one of the options ("saliva"; "nasal", or "others") for the type of test performed using the self-test kits.
 - Step 6: On the next page, pick one of the choices for submitting the self-test results ("positive"; "negative", or "invalid").

DIRECTIONS FOR USE

For video demonstration, please scan the QR code below:


1. Preparation Prior to Testing

- Choose a place where the test can be conducted without interruption for 15-30 minutes. Equilibrate the test midstream by bringing it to room temperature {15-30°C (59F-86F)} for 15-30 minutes prior to testing.
- Wash your hands with soap and water for at least 20 seconds before testing. If soap and water not available, use hand sanitizer with at least 60% alcohol.
- DO NOT eat, drink, smoke, or chew gum for 30 minutes before saliva collection procedure.
- Severe mouth ulcers and bronchitis may affect the saliva collection. Avoid from contact the infection site from different test ones.
- Use running water to clean mouth 30 minutes before saliva collection.
- Samples that are heavily contaminated by oral food residues cannot be used.
- Open your test kit, and you should have:

- Take out the midstream.


2. Sample Application

- Insert the absorbent tip into the mouth. Make sure midstream is horizontally placed.
- Swab the absorbent tip in the mouth and tongue to collect saliva.
- Note:
 - During saliva collection, gently hold it in mouth and let saliva naturally adsorb on the absorbent tip.
 - DO NOT eat, drink, or smoke prior to testing for at least 30 minutes.
 - Any saliva sample is appropriate for testing, but saliva sample collected in the morning, before rinsing mouth, eating, or drinking, is recommended.


3. Application Time

- Take the absorbent tip out from the mouth when the saliva sample is seen moving across the result window in the center of the midstream.

- Wait for 10 minutes and read the results.

- Note:
 - The length of the application time is determined by the amount of saliva of the individual. The more the saliva, the shorter the time, and the less the saliva, the longer it takes.

4. Disposal

- Discard the test midstream into the biohazard waste bag and dispose it.

INTERPRETATION OF RESULTS

NEGATIVE RESULT:

One colored line appears in the control line region (C). No line appears in the test line region (T). A negative result indicates that SARS-CoV-2 antigen is not present in the sample or present below the detectable level of the test.

POSITIVE RESULT:

Two lines appear. One colored line should be in the control line region (C) and another apparent colored line should be in the test line region (T). A positive result indicates that SARS-CoV-2 antigen is detected in the sample.

- If positive result is obtained, notify Ministry of Health by either using self-notification at MySejahtera mobile application or notify the District Health Office.
- Covid-19 positive individual should contact the nearest COVID-19 Assessment Center (CAC) (<http://covid-19.moh.gov.my/hotline>) to inform and get advice for further action to be taken.

INVALID RESULT:

Control line region (C) fails to appear. Insufficient sample volume or incorrect procedural techniques are most likely reasons for control line region (C) failure. Review the procedure and repeat the test with a new test kit. If problem persists, discontinue using the test kit immediately and contact your local distributor.

NOTE:

The intensity of the color in test line region (T) will vary depending on the concentration of SARS-CoV-2 antigen present in the sample. Therefore, any shade of color in the test line region (T) should be considered positive.

LIMITATIONS

- The test result obtained from SARS-CoV-2 Antigen Saliva Lolly Test cannot be used as a sole indicator of clinical indication. The infection should be confirmed by a specialist along with other laboratory results, clinical symptoms epidemiology, and additional clinical data.
- Tests are presumptive only and any positive results need to be confirmed by a laboratory PCR test and for follow-up clinical care.
- A negative test result does not rule out other viral or bacterial infections.
- If symptomatic and a negative result is obtained this should be confirmed immediately by laboratory PCR test.
- A false negative test can result if the amount of antigen in a sample is below the detection limit of the test or if the sample was collected incorrectly.
- Repeat antigen rapid testing is recommended every 24 hours for 3 days if there is a suspicion of infection, exposure to high-risk settings or other occupational risks.
- Excess blood or mucus on the specimen may interfere with test performance and may yield a false positive result.
- A positive test result for COVID-19 does not preclude an underlying co-infection with another pathogen, therefore the possibility of an underlying bacterial infection should be considered.
- Results from antigen testing should not be used as the sole basis to diagnose or exclude SARS-CoV-2 infection or to inform infection status.
- A positive test result does not rule out other viral or bacterial infections.
- The SARS-CoV-2 Antigen Saliva Lolly Test was evaluated on a small number of clinical samples. However, clinical performance is expected to be indicative of the predominant variations in circulation at the time and place of the clinical evaluation. Performance may vary based on the circulating variations, including novel SARS-CoV-2 strains and their predominance.

PERFORMANCE CHARACTERISTICS

Clinical Verification

The clinical performance of the SARS-CoV-2 Antigen Saliva Lolly Test was established with 232 sample collected from symptomatic patients, who with symptoms onset within 7 days.

| SARS-CoV-2 Antigen Saliva Lolly Test | Comparative RT-PCR Test Result | | |
|--------------------------------------|--------------------------------|--------------|-------|
| | Positive (+) | Negative (+) | Total |
| Detected Positive | 108 | 1 | 109 |
| Detected Negative | 7 | 116 | 123 |
| Total | 115 | 117 | 232 |
| Sensitivity | 93.91%, 95% CI (87.97,97.02) | | |
| Specificity | 99.15%, 95% CI (95.32, 99.85) | | |
| Accuracy | 96.55%, 95% CI (93.34, 98.24) | | |

Positive results broken down by days since symptoms onset.

| Days Since Symptoms Onset | RT-PCR Positive (+) | SARS-CoV-2 Antigen Saliva Lolly Test | PPA |
|---------------------------|---------------------|--------------------------------------|--------|
| 1 | 13 | 13 | 100% |
| 2 | 32 | 32 | 100% |
| 3 | 52 | 51 | 98.08% |
| 4 | 69 | 67 | 97.10% |
| 5 | 86 | 83 | 96.51% |
| 6 | 102 | 97 | 96.00% |
| 7 | 115 | 108 | 93.91% |

Positive results broken down by Ct values.

| SARS-CoV-2 Antigen Saliva Lolly Test | Comparative RT-PCR Test Result (Positive by Ct Value) | |
|--------------------------------------|---|------------------|
| | Positive (Ct≤25) | Positive (Ct<25) |
| Detected Positive | 69 | 39 |
| Total | 70 | 45 |
| Positive Agreement | 98.57% | 86.67% |

Detection Limit

When the virus is greater than 400TCID₅₀/mL, the positive detection rate is greater than 95%. When the virus content is less than 200TCID₅₀/mL, the positive detection rate is less than 95%. Thus, the minimum detection limit of SARS-CoV-2 Antigen Saliva Lolly Test is 400TCID₅₀/mL.

Precision

Three consecutive bathes of SARS-CoV-2 Antigen Saliva Lolly Test were tested for precision. Different batches of SARS-CoV-2 Antigen Saliva Lolly Test were used to test the same negative sample 10 times in succession, and the results were all negative. Different batches of SARS-CoV-2 Antigen Saliva Lolly Test were used to test the same positive samples 10 times in succession, and the results were all positive.

HOOK Effect

When the virus content in the sample to be tested reaches 4.0 x 10⁵ CID₅₀/mL, the test result still does not show the HOOK effect.

Cross-Reactivity

Cross-reactivity of the SARS-CoV-2 Antigen Saliva Lolly Test was evaluated. The results showed no cross-reactivity with the following sample.

| Name | Concentration |
|------------------------------|--|
| HCOV-HKU1 | 10 ⁵ TCID ₅₀ /ml |
| Staphylococcus aureus | 10 ⁶ TCID ₅₀ /ml |
| Group A streptococci | 10 ⁶ TCID ₅₀ /ml |
| Measles virus | 10 ⁵ TCID ₅₀ /ml |
| Mumps virus | 10 ⁵ TCID ₅₀ /ml |
| Adenovirus type 3 | 10 ⁵ TCID ₅₀ /ml |
| Mycoplasmal pneumonia | 10 ⁶ TCID ₅₀ /ml |
| Parainfluenzavirus.type2 | 10 ⁵ TCID ₅₀ /ml |
| Human metapneumovirus | 10 ⁵ TCID ₅₀ /ml |
| Human coronavirus OC43 | 10 ⁵ TCID ₅₀ /ml |
| Human coronavirus 229E | 10 ⁵ TCID ₅₀ /ml |
| Bordetella parapertusis | 10 ⁵ TCID ₅₀ /ml |
| Influenza B Victoria STRAIN | 10 ⁵ TCID ₅₀ /ml |
| Influenza B YSTRAIN | 10 ⁵ TCID ₅₀ /ml |
| Influenza A H1N1 2009 | 10 ⁵ TCID ₅₀ /ml |
| Influenza A H3N2 | 10 ⁵ TCID ₅₀ /ml |
| H7N9 | 10 ⁵ TCID ₅₀ /ml |
| H5N1 | 10 ⁵ TCID ₅₀ /ml |
| Epstein-Barr virus | 10 ⁵ TCID ₅₀ /ml |
| Enterovirus CA16 | 10 ⁵ TCID ₅₀ /ml |
| Rhinovirus | 10 ⁵ TCID ₅₀ /ml |
| Respiratory syncytial virus | 10 ⁵ TCID ₅₀ /ml |
| Streptococcus pneumoni-ae | 10 ⁶ TCID ₅₀ /ml |
| Candida albicans | 10 ⁶ TCID ₅₀ /ml |
| Chlamydia pneumoniae | 10 ⁶ TCID ₅₀ /ml |
| Bordetella pertussis | 10 ⁶ TCID ₅₀ /ml |
| Pneumocystis jiroveci | 10 ⁶ TCID ₅₀ /ml |
| Mycobacterium tubercu- losis | 10 ⁶ TCID ₅₀ /ml |
| Legionella pneumophila | 10 ⁶ TCID ₅₀ /ml |
| Human coronavirus NL63 | 10 ⁵ TCID ₅₀ /ml |
| MERS coronavirus | 10 ⁵ TCID ₅₀ /ml |

Interfering Substances

The test results do not interfered with the substances at the following concentration.

| Interfering Substances | Conc. | Interfering Substances | Conc. |
|------------------------|--------|---------------------------------------|----------|
| Whole Blood | 4% | Compound Benzoin Gel | 1.5mg/ml |
| Ibuprofen | 1mg/ml | Cromolyn glycate | 15% |
| Tetracycline | 3ug/ml | Chloramphenicol | 3ug/ml |
| Mucin | 0.5% | Mupirocin | 10mg/ml |
| Erythromycin | 3ug/ml | Oseltamivir | 5mg/ml |
| Tobramycin | 5% | Naphazoline Hydrochloride Nasal Drops | 15% |
| Menthol | 15% | Fluticasone propionate spray | 15% |
| Afrin | 15% | Deoxyepinephrine hydrochloride | 15% |

FREQUENTLY ASKED QUESTIONS

1. Will this test hurt?

No, the flat pad is not sharp and it should not hurt. Sometimes the tongue can feel slightly uncomfortable or dry. If you feel pain, please stop the test and seek advice from a healthcare provider.

2. How do I know that the test was run properly?

A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms adequate membrane wicking.

3. What should I do if the result shows positive?

You must take a laboratory PCR test immediately, self-isolate and contact your local health authority for further advice. You should also inform the immediate contacts you have had in past 24 hours so they can take any appropriate precautions.

4. What should I do if the result shows negative?

Negative results may require additional testing to confirm your results if you are symptomatic. If symptomatic, continue antigen testing every 24 hours for 3 days or take a laboratory PCR test. If asymptomatic, it is likely that you were not infectious at the time the test was taken. A negative test result, however, is not a guarantee that you do not have coronavirus. Please continue to follow social distancing, washing hands regularly and wearing masks as directed.



5. Can Sejoy SARS-CoV-2 Antigen Saliva Lolly Test detect various variants of COVID-19?

Yes, Sejoy SARS-CoV-2 Antigen Saliva Lolly Test can detect Alpha, Beta, Gamma, Delta and Omicron COVID-19 mutants based on the studies conducted so far.

BIBLIOGRAPHY

- Weiss SR,Leibowitz JZ.Coronavirus pathogenesis. Adv Virus Res 2011;81:85-164
- Cui J,Li F,Shi ZL.Origin and evolution of pathogenic coronaviruses.Nat Rev Microb iol 2019;17:181-192.
- Su S,Wong G,Shi W,et al.Epidemiology.genetic recombination,and pathogenesis of coronaviruses. TrendsMicrobiol 2016;24:490-502.

Index of Symbols

| | | | | | |
|---|----------------------------------|---|---------------|---|----------------------------------|
|  | Consult Instruction for use |  | Tests per kit |  | Do not use if package is damaged |
|  | For in vitro diagnostic use only |  | Use by date |  | Do not reuse |
|  | Store between 2-30°C |  | Lot Number |  | Catalogue number |
|  | Keep away from sunlight |  | Keep dry | | |

Hangzhou Sejoy Electronics & Instruments Co.,Ltd.

Area C, Building 2, No.365, Wuzhou Road, Yuhang Economic Development Zone, 311100 Hangzhou City, Zhejiang, China
Website: www.sejoy.com

Authorized Representative (Malaysia):

Medical Innovation Ventures Sdn. Bhd. (MEDIVEN)
1st Floor, Plot 88f, Lintang Bayan Lepas 10, Bayan Lepas Industrial Park, Phase 4, 11900 Bayan Lepas, Pulau Pinang, Malaysia.
Tel: +604 305 2730 Fax: +604 305 2730
Website: www.mediven.com.my
Email: sales@mediven.com.co

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Hangzhou Sejoy Electronics & Instruments Co., Ltd.